


5'-Aminoxy-modifier-11 CEP

Product No. BA 0350

Product Information



$C_{38}H_{54}N_3O_8P$
Mol. Wt.: 711.82

The use of oxime formation in ligation reactions of oligonucleotides has been widely established in the literature for over a decade.¹ As examples, 5'-aminoxy-modifiers have been used in oxime ligation for peptide-oligonucleotide conjugates,² attachment of nucleosides to solid supports,³ and head to tail cyclization of oligonucleotides.⁴ For the synthesis of an oligo bearing a 5'-aminoxy group we now offer 5'-Aminoxy-modifier-11 CEP (**BA 0350**).

Use: For oligonucleotide synthesis, employ acetonitrile diluent at the concentration recommended by the synthesizer manufacturer. Use standard coupling protocols; in our hands, extended coupling times were not required, and coupling efficiencies of > 99 % could be obtained. Cleavage from the solid support may be carried out by standard procedures. Standard nucleobase deprotection conditions may be employed.

References

1. For recent reviews see: a) Zatsepin, T.S.; Stetsenko, D.A.; Gait, M.J.; Oretskaya, T.S. *Bioconjugate Chem.*, **2005**, *16*(3), 471-89. b) Singh, Y.; Edupuganti, O.P.; Villen, M.; Defrancq, E.; Dumy, P. *Comptes Rendus Chim.*, **2005**, *8*(5), 789-96.
2. a) Cebon, B.; Lambert, J.N.; Leung, D.; Mackie, K.; McCluskey, K.L.; Nguyen, H.; Tassone, C. *Aust J. Chem.*, **2000**, *53*, 3333-40. b) Prater, C.E.; Miller, P.S. *Bioconjugate Chem.*, **2003**, *14*(2), 320-30. c) Prater, C.E.; Miller, P.S. *Bioconjugate Chem.*, **2004**, *15*(2), 498-507.
3. a) Salo, H.; Virta, P.; Hakala, H.; Prakash, T.P.; Kawasaki, A.M.; Manoharan, M.; Lönnberg, H. *Bioconjugate Chem.*, **1999**, *10*(5), 815-23. b) Defrancq, E.; Hoang, a.; Vinet, F.; Dumy, P.; *Bioorg. Med. Chem. Lett.*, **2003**, *13*, 2683-6. c) Bincheva, M.; Scheibler, L.; Lincoln, P.; Vogel, H.; Akerman, B. *Langmuir*, **1999**, *15*, 4317-20.
4. Edupuganti, O.P.; Defrancq, E.; Dumy, P. *J. Org. Chem.*, **2003**, *68*, 8708-10.