5-Aminoallyl-dU CEP Product No. BA 0311 Product Information



Allows the introduction of a 5-aminoallyl-dU residue into oligonucleotides for the purpose of post-synthetic labeling by acylation. Further, triplex-forming oligonucleotides (TFOs) containing this residue form triplexes that are similar in stability to those bearing unmodified residues.

Early work on the introduction of a 5-aminoallyl-dU amino modifier into oligonucleotides involved the methyl phosphoramidite.¹ 5-Aminoallyl-dU CEP, employing a 2-cyanoethyl phosphoramidite, was reported later.^{2,3}

Williams and co-workers³ studied the incorporation of 5-aminoallyl-dU residues into triplex-forming oligonucleotides (TFOs), where a single incorporation did not change the stability of the triplex vs. unmodified TFOs. Multiple incorporations led to a lower triplex stability. The pKa of the protonated form of the amino group of 5-aminoallyl-dU is expected to be 9.7 based on studies on the free nucleoside. Conformational studies on the nucleoside were also carried out.

Note: We also offer 5-Aminoallyl-U CEP (BA 0269), the ribose version of this amine modifier, .

Use: Lermer, *et al.*,² suggest a coupling time of 10 min. Oligonucleotides containing 5aminoallyl-dU residues were deprotected with concentrated ammonium hydroxide at 55 °C overnight. Williams and co-workers³ used a phosphoramidite concentration of 1.5x normal. They also employed N^4 -acetyl-dC phosphoramidite followed by deprotection with concentrated ammonium hydroxide at room temperature for 12 h.

In our hands, acetonitrile diluent at the concentration recommended by the synthesizer manufacturer was employed. A standard coupling time was sufficient to obtain high yields; extending to 15 minutes led to only a marginal improvement. Cleavage from the solid support was carried out using standard procedures. Deprotection was carried out using concentrated ammonium hydroxide at 55 °C for 16-18 h.

BERRY&ASSOCIATES

Literature:

- 1. Cook, A. F.; Vuocolo, E.; Brakel, C. L. Nucleic Acids Res. 1988, 16, 4077-4095.
- 2. Lermer, L.; Yoann, R.; Ting, R.; Perrin, D. M. J. Am. Chem. Soc. 2002, 124, 9960-9961. See especially the Supporting Information.
- 3. Brazier, J. A.; Shibata, T.; Townsley, J.; Taylor, B. F.; Frary, E.; Williams, N. H.; Williams, D. M. *Nucl. Acids Res.* **2005**, *33*, 1362-1371.