2-Thiouridine CEP Product No. BA 0415

Product Information



C₄₅H₆₁N₄O₈PSSi Mol. Wt.: 877.12

The effect of 2-thiouridine on sugar conformation and RNA duplex thermodynamics has been well documented.^{1,2} The presence of the 2-sulfur modification stabilizes the 3'-endo sugar conformation at the nucleoside and nucleotide level.¹ Wobble base pair specificity can also be improved by substituting 2-thiouridine for uridine. Testa and co-workers² have shown the at S²U favors S²U-A pairing more than S²U-G pairing, and more than U favors U-A relative to U-G. The sulfur modification improves specificity while retaining other key uridine activities.

In addition, thiolated uridine has been shown to improve the rate and fidelity of both nonenzymatic³ and ribozyme⁴ catalyzed nucleotide addition in RNA synthesis. Our 2-Thiousidine CEP (BA 0415) allows for the efficient insertion of 2-thiouridine into oligonucleotides.

Use: Dissolve the phosphoramidite in acetonitrile at concentrations recommended by the synthesizer manufacturer. Coupling should be carried out using standard instrument RNA protocols. Improved yields may be obtained with a modified oxidation protocol.^{1,5} Cleavage from the solid support can be carried out under standard conditions, and best results are obtained when deprotection is done overnight at room temperature with ammonium hydroxide.

BERRY&ASSOCIATES

Rev. 17-April-2017 Vers. 1.0

References:

- 1. Kumar, R.K.; Davis, D.R. Nuc. Acid. Res. 1997, 25(6), 1272-1280.
- 2. Testa, S.M.; Disney, M.D.; Turner, D.H.; Kierzek, R. *Biochem.* **1999**, *38*, 16655-16662.
- 3. Heuberger, B.D.; Pal, A.; Del Frate, F.; Topkar, V.V.; Szostak, J.W. *J. Am. Chem. Soc.*, **2015**, *137* (7), pp 2769–2775.
- 4. Prywes, N.; Michaels, Y.S.; Pal, A.; Oh, S.S., Szostak, J.W. *Chem. Commun.* **2016**. *52*, 6529-6532.
- a) Kumar, R.K.; Davis, D.R. J. Org. Chem. 1995, 60, 7726-7727. b) Gait, M.J. Oligonucleotide Synthesis-A Practical Approach 1990, Oxford University Press, London.