
# 2'-Deoxyxanthosine CEP (dX CEP) Product No. BA 0313 Product Information




 $\begin{array}{c} C_{56}H_{61}N_8O_{12}P\\ Mol. \ Wt.: \ 1069.10 \end{array}$ 

2'-Deoxyxanthosine CEP (dX CEP) is useful for the installation of 2'deoxyxanthosine (dX) nucleotides into oligonucleotides.<sup>1,2</sup>

**Introduction:** Xanthosine (X) and 2'-deoxyxanthosine (dX) have been incorporated into oligoribonucleotides<sup>3</sup> and oligodeoxyribonucleotides,<sup>1,2,4-8</sup> imparting a variety of properties to these nucleic acids. Examples include:

• The pKa of the proton at  $N^3$  of dX is 5.5-5.7,<sup>2,8</sup> therefore dX is present as a mono-anion to the extent of about 95% at neutral pH.



• dX residues have the same hydrogen-bond donor/acceptor pattern as thymidine, thus allowing base-pairing with dA.<sup>8</sup> This leads to a distortion of the helix, since two purines are paired. Stronger base stacking interactions are apparent. Duplex stability variations as a function of the opposing base pair have been studied.<sup>4,8</sup>

## **BERRY**&ASSOCIATES

01-28-2008 Vers. 1.0

- dX residues are a product of deamination of dG.<sup>2,8,9</sup>
- dX is susceptible to depurination under acidic conditions, forming abasic sites.<sup>2,8</sup> At pH 4, dX depurinates more quickly than dG by a factor of about 1000. At neutral pH, dX has about the same stability as dG. The recognition and excision of dX by base excision repair enzymes has been studied.<sup>2</sup>
- dX has been considered to be a promiscuous or universal nucleoside.<sup>1,8</sup>
- dX may expand the genetic code when paired to a pyrimidine-2,4diamine C-nucleoside.<sup>6</sup>
- Studies of polymerase-mediated incorporation of various nucleotides opposite dX have been carried out.<sup>2,4</sup>
- The presence of dX residues may lead to triplexes, quartets, and other G-stranded structures.<sup>8</sup>

For the incorporation of dX into oligonucleotides, we have chosen 2'-Deoxyxanthosine CEP (BA 0313), which bears 2-(4-nitrophenyl)ethyl (NPE) protecting groups on  $O^2$  and  $O^{6,1,2}$  Other monomers for the incorporation of dX or X have been reported.<sup>1,3,4-8</sup> Those without  $O^2/O^6$  protection appear to lead to byproducts in their synthesis and use. A recent alternative to NPE protection involves allyl groups, which may be removed via Pd(0) catalysis.<sup>8</sup>

# Coupling, deprotection, and purification:

Benner, Battersby, and co-workers<sup>1</sup> synthesized oligonucleotides using dX CEP and standard methods with the exception of an extended coupling time of 3 minutes, the use of 1*H*-imidazole-4,5-dicarbonitrile (0.25 M) as activator, and dichloroacetic acid (2.5% in dichloromethane) to detritylate. The alternate activator and deblock reagent were precautions to avoid depurination of dX. For cleavage and deprotection, the solid support bearing the trityl-off oligonucleotide was treated with concentrated aqueous ammonium hydroxide for 24-40 h at rt. Evaporation, suspension of the residue in water, filtration through a 0.2 µm filter, and evaporation was carried out. The residue was dissolved in and aqueous dioxane solution containing *N*,*N*,*N*,*N*-tetramethylguanidinium 2nitrobenzaldoximate and heated at 70 °C for 18 h, followed by quenching with sodium citrate and various extractions (see reference 1 for details). If the cleavage and deprotection was carried out with ammonium hydroxide only (48 h at 60 °C), about a third of the oligonucleotide had iso-dG residues, resulting from displacement of the 2-(4-nitrophenyl)ethoxy group at position 6 of the purine nucleobase.

*Wuenschell, O'Connor, and Termini*<sup>2</sup> synthesized oligonucleotides using dX CEP via standard methods, except for the presence of 4 Å molecular sieves in the acetonitrile solution of the phosphoramidite. The oligonucleotide was synthesized trityl-off, then cleaved from the support with 30% aqueous ammonia at 50 °C for 16 h. The NPE protecting groups were subsequently removed with DBU in dry pyridine at rt for 18 h. Various additional manipulations (see reference 2 for details) gave the desired material. No mention of iso-dG formation was made.

In our hands, dissolving 2'-Deoxyxanthosine CEP in dry acetonitrile without sieves followed by standard protocols as recommended by the synthesizer manufacturer led to >95% yield of incorporation of dX nucleotides with no requirement for extended coupling times. The only modification was the use of 0.25 M 5-ethylthio-1*H*-tetrazole (ETT) rather than 1*H*-tetrazole, the use of which was not attempted. Cleavage from the controlled-pore glass solid support was carried out with concentrated aqueous ammonium hydroxide at rt as per normal protocols. The resultant ammonium hydroxide solution can be used to deprotect other (non-dX) nucleobase protecting groups (e.g., rt for 24 h or 55 °C for 16 h; see Benner, *et al.*, for concerns about higher temperatures<sup>1</sup>). Under these conditions, we found that the NPE groups were not removed. Consult references 1 and 2 for protocols for NPE removal. Finally, note that the two NPE groups impart substantial lipophilicity to the oligonucleotide, resulting in a longer retention time on an RP-HPLC column than normal.

#### **References:**

1. Jurczyk, S. C.; Horlacher, J.; Devined, K. G.; Benner, S. A.; Battersby, T. R. *Helv. Chim. Acta* **2000**, *83*, 1517-1524.

2. Weunschell, G. E.; O'Connor, T. R.; Termini, J. *Biochemistry*, **2003**, *42*, 3608-3616.

3. Tuschl, T.; Ng, M. M. P.; Pieken, W.; Benseler, F.; Eckstein, F. *Biochemistry* **1993**, *32*, 11658-11668.

4. Eritja, R.; Horowitz, D. M.; Walker, P. A.; Ziehler-Martin, J. P.; Boosalis, M. S.; Goodman, M. F.; Itakura, K.; Kaplan, B. E. *Nucleic Acids Res.* **1986**, *14*, 8135-8153.

#### **BERRY**&ASSOCIATES

5. Van Aerschot, A.; Mag, M.; Herdewijn, P.; Vanderhaeghe, H. *Nucleosides Nucleotides Nucleic Acids* **1989**, *8*, 159-178.

6. Lutz, M. J.; Held, H. A.; Hottiger, M.; Hubscher, U.; Benner, S. A. *Nucleic Acids Res.* **1996**, *24*, 1308-1313.

7. Groebke, K.; Hunziker, J.; Fraser, W.; Peng, L.; Diederichsen, U.; Zimmermann, K.; Holzner, A.; Leumann, C.; Eschenmoser, A. *Helv. Chim. Acta* **1998**, *81*, 375.

8. Seela, F.; Shaikh, K. I. Helv. Chim. Acta 2006, 89, 2794-2814.

9. Suzuki, T.; Yamaoka, R.; Nishi, M.; Ide, H.; Makino, H. J. Am. Chem. Soc. 1996, 118, 2515.

## **BERRY**&ASSOCIATES