



6-Aza-dU CEP
Product No. BA 0303
Product Information



$C_{46}H_{52}N_5O_{10}P$
Mol. Wt.: 865.91

As compared to dU residues, 6-aza-dU nucleotides are deprotonated at neutral pH, have enhanced 3'-exonuclease stability, and exhibit a high-anti glycosidic bond conformation. Metal-DNA (M-DNA) complexes are formed at neutral pH.

6-Aza-2'-deoxyuridine (6-aza-dU)¹ has a nitrogen atom in place of the methine group at position 6 of the uracil ring. This results in a significant lowering of the pK_a of the N³ hydrogen.² Unlike dU, 6-aza-dU is deprotonated at neutral pH, thus bestowing the uracil base with a negative charge. Further, there are significant conformational changes imparted by the presence of the 6-aza substitution: 6-Aza-dU has a high-anti conformation and is present in solution with *N*-type sugar pucker.



Seela and Chittepu have reported that 6-aza-dU nucleotides may be incorporated into oligodeoxyribonucleotides using phosphoramidite chemistry.^{2,3} The N³-*o*-anisoyl protecting group⁴ is necessary to afford stability. DNA duplexes containing 6-aza-dU•dA base pairs have been studied² and are less stable than T•dA base pairs at neutral pH due to deprotonation of the N³ hydrogen on 6-aza-dU. However, at lower pH, duplex stability increases as N³ becomes protonated and therefore able to hydrogen bond to dA. The fact that 6-aza-dU residues are deprotonated at neutral pH results in easy metal-DNA (M-DNA) formation.² Canonical DNA forms M-DNA only at high pH, whereas 6-aza-dU-

containing DNA can form M-DNA at neutral or high pH. Finally, oligonucleotides containing 6-aza-dU show enhanced resistance to 3'-exonucleases (e.g., snake venom phosphodiesterase).

Coupling, deprotection, and purification: 6-Aza-dU CEP couples with greater than $\geq 95\%$ efficiency using the standard protocols recommended for popular synthesizers. Seela² recommends deprotection of oligomers with 25% aqueous ammonia at 60 °C for 16 h.⁵ Purification of the trityl-on oligonucleotide was accomplished by RP-HPLC (RP-18 column), after which the trityl group was removed with 2.5 % dichloroacetic acid in dichloromethane for 5 min. The detritylated oligomers were further purified via RP-HPLC and desalted on a short RP-18 silica column.

References:

1. (a) Kara, J.; Sorm, F. *Coll. Czech. Chem. Commun.* **1963**, 28, 1441-1448. (b) Holy, A.; Cech, D. *Coll. Czech. Chem. Commun.* **1974**, 39, 3157-3167. (c) Drasar, P.; Hein, L.; Beranek, J. *Coll. Czech. Chem. Commun.* **1976**, 41, 2110-2123. (d) Freskos, J. N. *Nucleosides Nucleotides* **1989**, 8, 549-555.
2. Seela, F.; Chittepu, P. *J. Org. Chem.* **2007**, 72, 4358-4366.
3. For the incorporation of 6-azauridine into oligoribonucleotides, see: (a) Beigelman, L.; Karpeisky, A.; Usman, N. *Nucleosides Nucleotides* **1995**, 14, 895-899 and (b) Oyelere, A. K.; Strobel, S. A. *Nucleosides Nucleotides Nucleic Acids* **2001**, 20, 1851-1858.
4. Welch, C. J.; Chattopadhyaya, J. *Acta Chem. Scand. Ser. B* **1983**, 37, 147-150.
5. In our hands, the N³-*o*-anisoyl protecting group was removed in a few minutes with concentrated ammonium hydroxide. Presumably, other standard cleavage/deprotection protocols are viable.