

Manual

ClearColi BL21(DE3) Electrocompetent Cells

For Research Use Only. Not for use in diagnostic procedures.

IMPORTANT
-80 °C storage required
immediately upon receipt

Manual

ClearColi BL21(DE3) Electrocompetent Cells

Contents

<u>1. Product description</u>	3
<u>2. Product specifications</u>	4
<u>3. Product designations and kit components</u>	4
<u>4. Storage conditions</u>	4
<u>5. Introduction to ClearColi technology</u>	5
5.1. Growth/colony characteristics of ClearColi BL21(DE3) Cells.....	5
<u>6. Endotoxin and ClearColi Cells</u>	6
6.1. LAL assay testing results.....	6
6.2. Lack of endotoxin in ClearColi Cells for mammalian cell applications	6
6.3. Avoiding endotoxin contamination from other sources.....	7
6.4. Measuring endotoxin prior to downstream applications.....	7
6.5. Endotoxin removal prior to downstream applications.....	7
<u>7. Preparation for transformation</u>	7
<u>8. Transformation protocol</u>	8
<u>9. Protein induction protocol</u>	9
<u>10. Media recipes</u>	10
<u>11. Technical support and product guarantee</u>	10

Manual

ClearColi BL21(DE3) Electrocompetent Cells

1. Product description

ClearColi™ BL21(DE3) cells are the first commercially available competent cells with a modified lipopolysaccharide (LPS) (Lipid IV_A - see figure 1) that does not trigger the endotoxic response in human cells. ClearColi cells lack outer membrane agonists for hTLR4/MD-2 activation therefore activation of hTLR4/MD-2 signaling by ClearColi cells is several orders of magnitude lower as compared with *E. coli* wild-type cells. Heterologous proteins prepared from ClearColi cells are virtually free of endotoxic activity. After minimal purification from ClearColi cells, proteins or plasmids (which may contain Lipid IV_A) can be used in most applications without eliciting an endotoxic response in human cells (figure 2).




Figure 1. The LPS of normal *E. coli* cells compared to the genetically modified lipid IV_A from ClearColi cells. In ClearColi cells, the oligosaccharide chain is deleted and two of the six acyl chains are removed to disable the endotoxin signal.




Figure 2. Comparison of endotoxic response from protein derived from ClearColi BL21(DE3) and traditional BL21(DE3) competent cells. ApoA1 protein was expressed from a T7-promoter based plasmid in both cell types, followed by a simple Ni-column purification step and without any subsequent endotoxin removal steps. The purified proteins were then tested for Toll-like receptor (TLR) stimulation by assessing NF-κB activation in HEK293 cells expressing human TLR4. Protein derived from ClearColi BL21(DE3) cells demonstrated no activation at concentrations five orders of magnitude greater than the protein from traditional BL21(DE3) cells.

Manual

ClearColi BL21(DE3) Electrocompetent Cells

2. Product specifications

Format:	DUO: 50 µL per vial, enough for two 25 µL reactions
Transformation efficiency:	Electrocompetent cells: $\geq 1 \times 10^9$ cfu/µg pUC19*
Genotype	<i>F⁻ ompT hsdS_B (r_B: m_B) gal dcm lon λ(DE3 [lacI lacUV5-T7 gene 1 ind1 sam7 nin5]) msbA148 ΔgutQ ΔkdsD ΔlpxL ΔlpxM ΔpagP ΔlpxP ΔeptA</i>
Expiration dating	ClearColi competent cells are guaranteed for performance for 1 year after receipt of the product. Do not use this product past the expiration date. If you experience technical performance issues or have general questions, please contact our Technical Support Team: techsupport@lgcgroup.com

3. Product designations and kit components

Strain	Kit size	Catalogue number	Reagent description	Part number	Transformations	Storage
ClearColi BL21(DE3) Electrocompetent Cells (White cap)	12 rxns (DUO)	60810-1	ClearColi BL21(DE3) Electrocompetent Cells	F862108	(6 x 50 µL)	-80 °C
			pUC19 Transformation Control	F92078-1	(1 x 20 µL)	-20 °C to -80 °C
			Expression Recovery Medium**	F98405-1	(1 x 12 mL)	-20 °C to -80 °C
	24 rxns (DUO)	60810-2	ClearColi BL21(DE3) Electrocompetent Cells	F862108	(12 x 50 µL)	-80 °C
			pUC19 Transformation Control	F88912-1	(1 x 20 µL)	-20 °C to -80 °C
			Expression Recovery Medium**	F98405-1	(2 x 12 mL)	-20 °C to -80 °C

* As a control for transformation, ClearColi Electrocompetent Cells are provided with supercoiled pUC19 DNA at a concentration of 10 pg/µl. Use 1 µL for transformation.

** Additional Expression Recovery Medium (lactose minus) can be ordered separately as catalogue # 80030-1, 96 mL (8 x 12 mL).

4. Storage conditions

All ClearColi Competent Cells require storage at -80 °C.

The cells are shipped on dry ice in one container, along with Expression Recovery Medium and supercoiled control pUC19 DNA.



Manual

ClearColi BL21(DE3) Electrocompetent Cells

5. Introduction to ClearColi technology

In ClearColi cells, two of the secondary acyl chains of the normally hexa-acylated LPS have been deleted, eliminating a key determinant of endotoxicity in eukaryotic cells. The six acyl chains of the LPS are the trigger which is recognized by the Toll-like receptor 4 (TLR4) in complex with myeloid differentiation factor 2 (MD-2), causing activation of NF- κ B and production of proinflammatory cytokines. The deletion of the two secondary acyl chains results in lipid IV_A, which does not induce formation of the activated heterotetrameric TLR4/MD-2 complex and thus does not trigger the endotoxic response. In addition, the oligosaccharide chain is deleted, making it easier to remove the resulting lipid IV_A from any downstream product.

ClearColi competent cells have a genetically modified LPS that disables the trigger for endotoxic response in mammalian cells. This was accomplished by incorporating seven genetic deletions that modify LPS to lipid IV_A, while one additional compensating mutation (*msbA148*) enables the cells to maintain viability in the presence of the LPS precursor lipid IV_A. ClearColi BL21(DE3) cells are BL21(DE3) derived cells with several key mutations which result in significantly reduced endotoxicity. The specific set of mutations made to ClearColi cells is as follows:

msbA148 Δ gutQ Δ kdsD Δ lpxL Δ lpxM Δ pagP Δ lpxP Δ eptA

5.1 Growth/colony characteristics of ClearColi BL21(DE3) Cells

ClearColi BL21(DE3) cells grow at approximately 50% of the rate of normal BL21(DE3) cells. Users should expect to see very small colonies for the first 24 hours after plating transformants. LGC Biosearch Technologies recommends incubating plates for 32-40 hours before picking colonies for future experiments (see Transformation protocol and Protein induction protocol sections of this manual for more information). Longer growth times are necessary to reach desired cell densities prior to inducing protein expression. In addition, users may observe some variation in colony size when plating ClearColi cells. A small portion (<2%) of colonies may be larger than the general population. These larger colonies have similar protein expression levels and endotoxin levels as the average size colonies.

5.2 Important information about medium formulation

ClearColi cells have a specially engineered membrane composition and require a modified medium formulation compared to most E. coli strains. Consider the following:

- *ClearColi cells are osmosensitive and require 1% NaCl in their growth medium. We strongly recommend you use high salt **LB-Miller Medium** to achieve optimal growth. **LB-Miller** differs from **LB-Lennox**. See section 10 for the LB-Miller recipe.*
- ***We do not recommend growing ClearColi cells in LB-Lennox or Super broth medium. These media have resulted in slow and suboptimal cell growth.***
- ***Do not include Mg²⁺ and Ca²⁺ in your medium. They have been shown to inhibit growth of ClearColi cells.***
- *ClearColi cells tend to aggregate when grown in liquid culture. We suggest that you vortex the cell solution before measuring OD₆₀₀.*

Manual

ClearColi BL21(DE3) Electrocompetent Cells

6. Endotoxin and ClearColi Cells

6.1 LAL assay testing results

Limulus amoebocyte lysate (LAL) assay testing is an FDA-approved method for detection of endotoxins and is the most common assay used. However, the LAL assay is an inappropriate method to discriminate between endotoxically active hexa-acylated LPS and endotoxically inactive tetra-acylated lipid IV_A. The structural requirements for activation of the LAL cascade by endotoxins differ from those for stimulation of the human immune cell system. While the acylation pattern of LPS/lipid A is a key determinant for stimulation of human immune cells, activation of the LAL cascade only marginally, if at all, depends on the number of acyl chains. Instead, reactivity in the LAL assay requires the 4'-monophosphoryl-diglycosamine backbone structure, which is present in both hexa-acylated LPS and tetra-acylated lipid IV_A of *E. coli*. As such, false positive results are due to the lack of specificity of the assay. The LAL assay recognizes a wider spectrum of LPS/lipid A variants than the central cellular endotoxin sensor system of the human immune cell system.

A simple Ni-column purification step for proteins produced from ClearColi cells significantly reduces LAL response levels. For example, Biosearch Technologies demonstrated a 99% reduction in LAL response comparing ApoA1 expressed using ClearColi cells versus *E. coli* EXPRESS BL21(DE3) Electrocompetent Cells.

Cell line	LAL results (endotoxin units/mg)	Percent reduction
ClearColi Electrocompetent Cells	450	99.1%
<i>E. coli</i> EXPRESS BL21(DE3) Electrocompetent Cells	53,800	

Residual endotoxin unit (EU) measurements are due to the non-specific nature of the assay unless extraneous LPS contamination from other sources is present. Alternative endotoxicity assays, such as those using HEK-Blue cells (InvivoGen), suggest that even in the presence of EU levels above normal thresholds targeted by researchers the actual stimulating effects from ClearColi-derived proteins are non-existent. Due to the non-specific nature of the LAL assay giving false positive endotoxic results with lipid IV_A from ClearColi cells, it is suggested that researchers consider alternative physiologically relevant methods. These methods include measurement of TLR stimulation as assessed by NF- κ B activation in HEK293 cells or human macrophage assays.

6.2 Lack of endotoxin in ClearColi Cells for mammalian cell applications

Lipid IV_A from ClearColi cells is incapable of inducing an endotoxic response in human immune cells. Seven separate deletions ensure that ClearColi cells cannot revert to normal LPS production. With proper controls, plasmids and proteins can be produced from ClearColi cells without need for downstream endotoxin removal steps. However, LPS contamination may be prevalent in your laboratory, and care must be taken to minimize LPS sources other than your cell strain.

While lipid IV_A is known as an endotoxin antagonist in human LPS-responsive cells, it has to be taken into account that the tetra-acylated LPS precursor may act as an endotoxic activator in other mammalian hosts such as mouse, Chinese hamster or equine cells, which reflects animal species-specific recognition and stimulatory activity of lipid IV_A due to species-specific differences in the structures of TLR4 and MD-2.

Manual

ClearColi BL21(DE3) Electrocompetent Cells

6.3 Avoiding endotoxin contamination from other sources

Although ClearColi cells will not produce endotoxin, it is still possible to contaminate your end product with endotoxins from other sources. Good laboratory sterile technique can adequately control extraneous LPS contamination. Biosearch Technologies recommends the following precautions:

- Use disposable pipette tips and centrifuge tubes certified as sterile and non-pyrogenic.
- Depyrogenate any glassware by heat treating at >250 °C for 1 hour prior to use.
- Do not use purification columns or resins that have come in contact with *E. coli*.
- Use reagents certified as low endotoxin or test reagents prior to use.
- Use a water source that is regularly tested for endotoxin contamination.
- Clean all laboratory surfaces with disinfectants.

6.4 Measuring endotoxin prior to downstream applications

In applications where minimal endotoxin levels are critical, it is strongly recommended that all normal precautions are taken. **Biosearch Technologies cannot guarantee a total absence of LPS due to the possibility of contamination from other sources.** Safety and downstream applications are the sole responsibility of the user.

6.5 Endotoxin removal prior to downstream applications

The need for endotoxin removal steps will depend on the user's method of endotoxin measurement and application. As previously discussed, the lipid IV_A of ClearColi cells does not cause an endotoxic response in human cells; however, the use of LAL testing may result in a relatively low EU measurement. Normally, a simple plasmid purification or Ni-column protein purification will be sufficient to lower LAL levels below threshold. If lower levels are desired, additional cleanup steps should be taken.

7. Preparation for transformation

ClearColi Electrocompetent Cells are provided in 50 µL aliquots (DUOs), sufficient for two transformation reactions of 25 µL each. Transformation is carried out in a 0.1 cm gap cuvette. Optimal settings for electroporation are listed in the table below. Typical time constants are 3.5 to 4.5 msec.

Optimal setting	Alternate settings (~ 20-50% lower efficiencies)	For 2.0 mm cuvettes
1.0 mm cuvette 10 µF 600 ohms 1,800 volts	1.0 mm cuvette 25 µF 200 ohms 1,400-1,600 volts	2.0 mm cuvette 25 µF 750 ohms 2,400 volts

Suggested electroporation systems: Bio-Rad MicroPulser #165-2100; Bio-Rad *E. coli* Pulser #165-2102; Bio-Rad Gene Pulser II #165-2105; BTX ECM630 Electroporation System.

Suggested electroporation cuvettes: Successful results are obtained with cuvettes from BTX (Model 610), Bio-Rad (Cat. #165-2089), or Eppendorf (Cat. # 4307-000-569). Users have reported difficulties using Invitrogen cuvettes (Cat.# 65-0030)

Optional transformation control reactions include electroporation with 1 µL (10 pg) of supercoiled pUC19 DNA.

Manual

ClearColi BL21(DE3) Electrocompetent Cells

To ensure successful transformation results, the following precautions must be taken:

- For best results ligation reactions must be purified or heat killed at 70 °C for 15 minutes before transformation.
- The DNA sample to be used for electroporation must be dissolved in water or a buffer with low ionic strength, such as TE. The presence of salt in the electroporation sample leads to arcing at high voltage, resulting in loss of the cells and DNA.
***NOTE:** Ligation reactions performed with Biosearch Technologies' CloneDirect® Ligation Buffer (included with Biosearch Technologies' Cloning or Ligation Kits) can be used immediately after heat inactivation without purification of the ligation products.*
- Microcentrifuge tubes and electroporation cuvettes must be thoroughly pre-chilled on ice before use.
- The cells must be completely thawed **on ice** before use.
- For highest transformation efficiency, use the provided Expression Recovery Medium to resuspend the cells after electroporation. Use of TB or other media will result in lower transformation efficiencies.

8. Transformation protocol

Before proceeding with transformation, refer to section 5.2 regarding medium formulation.

Use the following protocol for cells provided in microcentrifuge tubes.

1. Have Expression Recovery Medium and 17 mm x 100 mm sterile culture tubes readily available at room temperature (one tube for each transformation reaction). Transformation efficiency may decrease with the use of SOC or other media.
2. Place electroporation cuvettes (0.1 cm gap) and microcentrifuge tubes on ice (one cuvette and one microfuge tube for each transformation reaction).
3. Remove ClearColi cells from the -80 °C freezer and place on wet ice until they thaw completely (10-20 minutes).
4. When cells are thawed, mix them by tapping gently. Aliquot 25 µL of ClearColi cells to the chilled microcentrifuge tube on ice.
5. Add 1 µL of DNA or heat-denatured ligation reaction to the 25 µL of cells on ice. (Failure to heat-inactivate the ligation reaction will prevent transformation.) Stir briefly with pipette tip; do not pipet up and down to mix, because this can introduce air bubbles and warm the cells. Use of more than 2 µL of ligation mix may cause electrical arcing during electroporation..
6. Carefully pipet 25 µL of the cell/DNA mixture into a chilled electroporation cuvette without introducing bubbles. Quickly flick the cuvette downward with your wrist to deposit the cells across the bottom of the well. Electroporate according to the conditions recommended above.
7. Within 10 seconds of the pulse, add 975 µL of Expression Recovery Medium to the cuvette and pipet up and down three times to resuspend the cells. Transfer the cells and Expression Recovery Medium to a culture tube.
8. Place the tube in a shaking incubator at 200-250 rpm for 1 hour at 37 °C.
9. Spread up to 100 µL of transformed cells on LB-Miller plates containing the appropriate antibiotic.
10. Incubate the plates 32-40 hours at 37 °C. Very small colonies may be visible at 24 hours.

Manual

ClearColi BL21(DE3) Electrocompetent Cells

NOTE: Do not store ClearColi strains on agar plates at 4 °C for more than one week. While cells remain viable, there is an increased lag time before growth. Transformed ClearColi strains can be stored at -80 °C in LB Miller medium with the addition of 20% glycerol.

11. Transformed clones can be further grown in in LB-Miller and the appropriate antibiotic (see section 10).

9. Protein induction protocol

1. Inoculate a single colony from a freshly streaked plate into 40 mL of LB-Miller medium containing the appropriate antibiotic for the plasmid and host strain. To minimise the amount of expression of the target protein prior to induction, add glucose to the growth medium at a concentration of 0.5% (w/v).
2. Incubate with shaking at 37 °C overnight.
3. Inoculate 1 L of LB-Miller medium (no glucose) containing the appropriate antibiotic with all 40 mL of the overnight culture prepared in step 2. Alternatively, measure the OD₆₀₀ of your overnight culture and inoculate to a final OD₆₀₀ of 0.1.
4. Incubate with continual shaking at 37 °C until the OD₆₀₀ reaches 0.6-0.8 (approximately 4-5 hours).
5. Add IPTG to a final concentration of 0.4-1 mM. (Prepare a 1 M solution of IPTG by dissolving 2.38 g of IPTG in water and adjust the final volume to 10 mL. Filter sterilize before use). To determine the optimal concentration of IPTG for maximum expression of the target protein, test a range of IPTG concentrations.
6. Incubate, shaking at 37 °C for 3-4 hours. To determine the optimal time for induction of the target protein, it is recommended that a time course experiment be performed varying the induction from 2-16 hours.

NOTE: Final OD₆₀₀ of the ClearColi cells may be as low as 50% of that normally achieved with standard BL21(DE3) cells due to slower growth rates.

7. Place the culture on ice for 10 minutes. Harvest cells by centrifugation at 5,000 x g for 10 minutes at 4 °C.
8. Remove the supernatant . Depending on your workflow and specific protein, you can proceed to:
 - a. Store the pellet at -20 °C.
 - b. Store the pellet at -80 °C.
 - c. Immediately continue with processing (isolation and purification of protein).

Manual

ClearColi BL21(DE3) Electrocompetent Cells

10. Media recipes

LB-Miller plates

Per liter: 5 g yeast extract
10 g tryptone
10 g NaCl
15 g agar

Add deionised water to 1 liter. Autoclave and cool to 55 °C before adding the appropriate filter-sterilised antibiotic (e.g. ≤ 30 mg/L kanamycin; ≤ 100 mg/L ampicillin or carbenicillin; ≤ 30 mg/L chloramphenicol).

Pour approximately 25 mL per Petri plate.

LB-Miller culture medium for growth of transformants

Per liter: 5 g yeast extract
10 g tryptone
10 g NaCl

Add all components to deionized water. Autoclave and cool to 55 °C.

11. Technical support and product guarantee

Biosearch Technologies is dedicated to the success and satisfaction of our customers. Our products are tested to assure they perform as specified when used according to our recommendations. It is imperative that the reagents supplied by the user are of the highest quality. Please follow the instructions carefully and contact our technical support representatives if additional information is necessary. We encourage you to contact us with your comments regarding the performance of our products in your applications. Thank you.

If you require any further support, please do not hesitate to contact our Technical Support Team: techsupport@lgcgroup.com.

Limited use license

If you are a commercial entity, your ClearColi kit purchase comes with a 12-month Limited Use License from Research Corporation Technologies.

Please refer to Legal Information (p. 11) or this URL for the terms of use: <http://clearcoli.com/licensing/limited-use-license-2/>

If you disagree with the terms of use, you have 10 days to contact Biosearch Technologies for permission to return unused cells.

Manual

ClearColi BL21(DE3) Electrocompetent Cells

Notice of limited label license, copyright, patents, warranties, disclaimers and trademarks

ClearColi Competent cells are subject to US Patent 8,303,964 and other US and foreign pending patents.

LGC Biosearch Technologies products are sold for research use only and are not to be used in humans or for medical diagnostics. Biosearch Technologies' liability with respect to any ClearColi product is limited to the replacement of the product. No other warranties of any kind, expressed or implied, including without limitation, any implied fitness for any particular use, are provided by Biosearch Technologies. Biosearch Technologies is not liable for any direct, indirect, incidental or consequential damages arising out of or in connection with the use or inability to use any of its ClearColi products.

Purchaser notification

Limited use label license

ClearColi System

The ClearColi System is a recombinant expression platform based on strains of *Escherichia coli* (*E. coli*) that are engineered to minimize human endotoxin signalling. Patents, strains (including derivatives and progeny) and other materials related to the ClearColi System and licenses for its use to express and/or clone products are owned by Research Corporation Technologies, Inc. ("RCT"), Tucson, Arizona. LGC Biosearch Technologies ("Biosearch Technologies") has a license from RCT to sell the ClearColi System and related information (collectively, "ClearColi Kits") to third-parties for non-commercial research purposes only. A separate license is required for any commercial use, including use of the ClearColi Kit for research purposes in product development, product production and/or product improvement. Before using the ClearColi Kit, please read the license agreement described below. If you do not agree to be bound by its terms, contact Biosearch Technologies within 10 days for authorization to return the unused ClearColi Kit to receive a full refund. If you do agree to the terms, please visit <http://clearcoli.com/registration/> and enter the requested registration information.

Biosearch Technologies grants you a non-exclusive license to use the enclosed ClearColi Kit for non-commercial research or evaluation purposes only. The ClearColi Kit is being transferred to you in furtherance of, and reliance on, such limited use license. You may not conduct any reverse engineering of the ClearColi System or use the ClearColi Kit for any Commercial Purpose (as defined below) without a license for such purpose from RCT, Tucson, Arizona. "Commercial Purpose" means: any use of the ClearColi Kit or ClearColi Products (as defined below) in a Commercial Product (as defined below); any use of the ClearColi Kit or ClearColi Products in the manufacture of a Commercial Product; any use of the ClearColi Kit or ClearColi Products in the performance of services for third parties; any sale or other commercialization of ClearColi Products; any use of the ClearColi Kit or ClearColi Products directly in the clinical development, improvement, or GMP manufacture of a Commercial Product; or the filing or preparation of any patent application claiming or describing, in whole or in part, the ClearColi System, any ClearColi Products, or any Commercial Products. "ClearColi Products" means products (including, but not limited to, proteins, peptides and DNA plasmids) expressed with the use of the ClearColi Kit, or with the use of any ClearColi strains (including derivatives and progeny). "Commercial Product" means any product intended for sale or commercial use. Commercial entities and non-profit entities engaged in commercial activities may conduct their research and evaluation for twelve (12) months from the original purchase shipping date at which time this license automatically terminates. Commercial entities and nonprofit entities engaged in commercial activities will be contacted by RCT during the evaluation period regarding their desire for a commercial license.

Access to the ClearColi Kit and ClearColi strains must be limited solely to those officers, employees and students of your entity who need access to perform the above-described research or evaluation. You must inform each such officer, employee and student of the provisions of this license agreement and require them to agree, in writing, to be bound by the provisions of this license agreement. You may not distribute the ClearColi Kit or any ClearColi strain contained herein or in the ClearColi Kit (or any derivative or progeny thereof), to others, even those within your own entity. You may only transfer modified, altered, or original ClearColi strains or material from the ClearColi Kit to a third party following written notification of, and written approval from, Biosearch Technologies and RCT so that the recipient can be properly licensed. You may not assign, sublicense, rent, lease or otherwise transfer this license agreement or any of the rights or obligations hereunder, except as expressly permitted by Biosearch Technologies and RCT.

This license agreement is effective until terminated. You may terminate it at any time by destroying all ClearColi strains (including derivatives and progeny) and ClearColi Products in your control. It will also terminate automatically if you fail to comply with the terms and conditions of the license agreement. You shall, upon termination of the license agreement, destroy all ClearColi strains and ClearColi Products in your control, and so notify Biosearch Technologies in writing. Information about commercial licenses may be obtained from Research Corporation Technologies, Inc. at the following address:

Brad Morie, MBA CLP,
Research Corporation Technologies, Inc.,
6440 N. Swan Rd, Suite 200,
Tucson, Arizona 85718.
Email: bmorie@rctech.com,
Tel: 520-748-4413, Fax: 520-748-0025.



X f in @LGCBiosearch | biosearchtech.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2024. All rights reserved. GEN/1013/EK/0924

For Research Use Only. Not for use in diagnostic procedures.



BIOSEARCH™
TECHNOLOGIES

GENOMIC ANALYSIS BY LGC