

Manual

RNase I

For Research Use Only. Not for use in diagnostic procedures.

RNase I is part of the Epicentre™ product line, known for its unique genomics kits, enzymes, and reagents which offer high quality and reliable performance.

Manual

RNase I

Contents

1. Introduction	3
2. Product designations and kit components	3
3. Product specifications	3
4. Protocol for removing RNA from DNA preparations	4
5. References	4
6. Further support	4

Manual

RNase I

1. Introduction

RNase I preferentially degrades single-stranded RNA to individual nucleoside 3' monophosphates by cleaving every phosphodiester bond.¹ By comparison, other ribonucleases cleave only after specific residues (e.g., RNase A cleaves 3' to pyrimidine residues). Thus, RNase I is useful for removing RNA from DNA preparations,² detecting mismatches in RNA:RNA and RNA:DNA hybrids^{2,3} and analysing and quantifying RNA in ribonuclease protection assays (RPA).^{4,5} The enzyme is completely inactivated by heating at 70 °C for 20 minutes in the presence of 5 mM dithiothreitol (DTT), eliminating the requirement to remove the enzyme prior to many subsequent procedures.

2. Product designations and kit components

Product	Kit size	Catalog number	Reagent description	Part number	Volume
RNase I, <i>E. coli</i>	1,000 Units	N6901K	RNase I (10 U/μL)	E0067-10D1	100 μL
			DTT (0.1 M)	SS000065-D1	2.5 mL
			RNase I Dilution Buffer	SS000255-D1	1 mL
			10X TNE Buffer	SS000806-D1	5 mL

3. Product specifications

Storage: Store only at -20 °C in a freezer without a defrost cycle.

Storage buffer: RNase I is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 0.1 M NaCl and 0.1 mM EDTA.

RNase I Dilution Buffer: A 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 0.1 M NaCl and 0.1 mM EDTA.

10X TNE Buffer: 100 mM Tris-HCl (pH 7.5), 1 M NaCl and 10 mM EDTA.

Unit definition: One unit degrades 100 ng of *E. coli* ribosomal RNA per second into acid-soluble nucleotides at 37 °C.

Quality control: RNase I is function-tested in a reaction containing 10 mM Tris-HCl (pH 7.5), 100 mM NaCl, 1 mM EDTA and 60 μg of *E. coli* ribosomal RNA with varying amounts of enzyme.⁶

Contaminating activity assays: RNase I is free of detectable exo- and endodeoxyribonuclease activities as judged by incubation of 1 μg of various DNA substrates with 4 x 10⁶ U of enzyme at 37 °C for 16 hours.

Manual

RNase I

4. Protocol for removing RNA from DNA preparations

RNase I can be used in place of RNase A for removing RNA from DNA preparations. In contrast to RNase A, RNase I effectively degrades contaminating RNA to mono- and dinucleotides that will not interfere with visualisation of small DNA molecules. After RNA removal, the enzyme can be inactivated by heating at 70 °C for 20 minutes in the presence of 5 mM DTT.

Protocol

1. Isolate DNA from 1-2 mL of overnight bacterial culture using a standard alkaline lysis procedure.⁵
2. After ethanol precipitation, suspend the DNA in 1X TNE buffer (page 3) at a concentration appropriate for subsequent applications (see Notes below).
3. Dilute RNase I enzyme 10-fold with RNase I Dilution Buffer and add 1.5-2 U to the DNA preparation.
4. Incubate at 37 °C for 30 minutes to degrade contaminating RNA.
5. Add DTT to a final concentration of 5-10 mM.
6. Incubate at 70 °C for 20 minutes to inactivate the enzyme.

Notes

Reaction buffer: Incubation with RNase I can be performed simultaneously with the digestion of plasmid DNA by restriction endonucleases. RNase I maintains ≥90% activity in buffers containing between 100 mM to 200 mM salt (either NaCl or KOAc). The activity of the enzyme is also relatively constant over a pH range of 7.0-8.8. Therefore, if the restriction endonuclease buffer is within these parameters, RNase I digestion can be performed in the restriction endonuclease buffer.

Enzyme dilution: Diluted enzyme may be stored for up to two months at -20 °C in a freezer without a defrost cycle.

5. References

1. Shen, V. and Schlessinger, D. (1982) *The Enzymes* **XV** (Part B), 501.
2. Winter, E. et al., (1985) *Proc. Natl. Acad. Sci. USA* **82**, 7575.
3. Myers, R.M. et al., (1985) *Science* **230**, 1242.
4. Sambrook, J. et al., (1989) in: *Molecular Cloning: A Laboratory Manual (2nd ed.)*, Cold Spring Harbor Laboratory Press, New York.
5. Saccomanno, C.F. et al., (1992) *BioTechniques* **13**, 847.
6. Corbishley, T.P. et al., (1984) *Meth. Enzymatic Anal.* **4**, 134.

6. Further support

If you require any further support, please do not hesitate to contact our Technical Support Team: techsupport@lgcgroup.com



Integrated tools. Accelerated science.



@LGCBiosearch

lucigen.com

biosearchtech.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2021. All rights reserved. GEN/889/EK/0121

For Research Use Only. Not for use in diagnostic procedures.



BIOSEARCH™ TECHNOLOGIES

GENOMIC ANALYSIS BY LGC