

NxGen[®] T4 DNA Ligase High Concentration Rapid Kit



IMPORTANT!

-20 °C Storage Required

Immediately Upon Receipt

FOR RESEARCH USE ONLY. NOT FOR HUMAN OR DIAGNOSTIC USE

Lucigen Corporation 2905 Parmenter St, Middleton, WI 53562 USA
Toll Free: (888) 575-9695 | (608) 831-9011 | FAX: (608) 831-9012
custserv@lucigen.com techserv@lucigen.com www.lucigen.com

NxGen® T4 DNA Ligase High Concentration Rapid Kit

Table of Contents

Technical Support	2
Product Description	2
Product Specifications	3
Product Designations and Kit Components	3
Components & Storage Conditions	3
Rapid Ligation Reaction Set-Up	4
Standard Ligation Reaction Set-Up	4
References	5
Notice of Limited Label License, Copyright, Patents, Warranties, Disclaimers and Trademarks	6

Technical Support

Lucigen is dedicated to the success and satisfaction of our customers. Our products are tested to assure they perform as specified when used according to our recommendations. It is imperative that the reagents supplied by the user are of the highest quality. Please follow the instructions carefully and contact our technical service representatives if additional information is necessary. We encourage you to contact us with your comments regarding the performance of our products in your applications. Thank you.

Lucigen Technical Support

Email: techserv@lucigen.com

Phone: (888) 575-9695

Product Guarantee: Lucigen guarantees that this product will perform as specified for one year from the date of shipment. Please avoid using reagents for greater than one year from receipt.

Product Description

T4 DNA Ligase catalyzes the formation of a phosphodiester bond between the terminal 5' phosphate and a 3' hydroxyl groups of duplex DNA or RNA. The enzyme efficiently joins blunt and cohesive ends and repairs single stranded nicks in duplex DNA, RNA, or DNA/RNA hybrids (1).

Storage buffer: T4 DNA Ligase is supplied in: 10 mM Tris-HCl, 50 mM KCl, 1 mM dithiothreitol, 0.1 mM EDTA, 0.1% Triton X-100, 50% glycerol, pH 7.5 @ 25 °C.

2X Rapid Ligation Buffer is composed of 132 mM Tris-HCl, 20 mM MgCl₂, 2 mM dithiothreitol, 2 mM ATP, 15% PEG, pH 7.6 @ 25 °C.

10X T4 DNA Ligase Buffer is composed of 500 mM Tris-HCl, 100 mM MgCl₂, 50 mM dithiothreitol, 10 mM ATP, pH 7.6 @ 25 °C.

Source of protein: A recombinant *E. coli* strain carrying the cloned T4 DNA Ligase gene.

NxGen® T4 DNA Ligase High Concentration Rapid Kit

Unit Definition: One Weiss unit is defined as the amount of enzyme required to convert 1 nmol of ³²P-labeled inorganic pyrophosphate into Norit absorbable material in 20 minutes at 37°C, using specified reaction conditions(2).

Note: 1 Weiss Unit is approximately 67 cohesive end units.

Product Specifications

TEST	SPECIFICATION
Purity (SDS-PAGE)	>99%
SS Exonuclease	6,000 U <0.1% released
DS Exonuclease	6,000 U <0.1% released
Endonuclease	6,000 U <0.1% converted
E. coli 16S rDNA Contamination	3,000 U <10 copies


Product Designations and Kit Components

Product	Ligase Concentration	Kit Size	Catalog Number	Reagent Description	Part Numbers	Volume
NxGen® T4 DNA Ligase (High Concentration)	10 Units* /μL	1,500 Units*	30243-1	T4 DNA Ligase (High Concentration)	F83914-1	150 μL
				10X T4 DNA Ligase Buffer	F88912-1	1.5 mL
				2X Rapid Ligation Buffer	F88915-1	1.5 mL
		7,500 Units*	30243-2	T4 DNA Ligase (High Concentration)	F83914-1	5 x 150 μL
				10X T4 DNA Ligase Buffer	F88912-1	5 x 1.5 mL
				2X Rapid Ligation Buffer	F88915-1	5 x 1.5 mL

*Weiss Units

Components & Storage Conditions

Store all Kits and Components at -20 °C



NxGen® T4 DNA Ligase High Concentration Rapid Kit

Rapid Ligation Reaction Set-Up

1. Add all of the components below to a clean reaction vessel.

Quantity	Component
10 µL	2X Rapid Ligase Buffer
X	Vector (1-10 ng/µL)
Y	Insert (1-10 ng/µL)
1 µL	T4 DNA Ligase (10 U/µL)
To 20 µL	Nuclease-free Water

2. Mix well by pipetting.
3. Incubate at 25 °C for 10 minutes
4. Purify DNA using a PCR clean-up column and elute in ~50 µL.

– OR –

Immediately dilute in TE or water (at least 1:10).

Note: Heat inactivation of ligation reactions containing 2X Rapid Ligation Buffer is not recommended, as it may reduce transformation efficiency.

5. Transform 0.1-10 ng of ligation product into a chemically or electrocompetent cell line that is compatible with the vector.

Notes:

- High-concentration T4 DNA Ligase in combination with the 2X Rapid Ligation buffer greatly stimulates the rate and efficiency of blunt-end ligation; therefore, incubations of over 10 minutes are NOT recommended and can greatly reduce the transformation efficiency. Use the recommended protocol to maximize transformation efficiency.
- 10X T4 DNA Ligase Buffer does not contain PEG and is compatible with Standard Ligation Reaction (below), includes longer incubation and heat inactivation.

Standard Ligation Reaction Set-Up

1. Add all of the components below to a clean reaction vessel.

Quantity	Component
2 µL	10X T4 Ligase Buffer
X	Vector (1-10 ng/µL)
Y	Insert (1-10 ng/µL)
1 µL	T4 DNA Ligase (10 U/µL)
To 20 µL	Nuclease-free Water

2. Mix well by pipetting.
3. Incubate at 25 °C for 30 minutes.
4. Heat inactivate the reaction by incubating the ligation at 70 °C for 15 minutes.

NxGen® T4 DNA Ligase High Concentration Rapid Kit

5. Purify DNA using a PCR clean-up column and elute in ~50 µL.

– OR –

Immediately dilute in TE or water (at least 1:10).

6. Transform 0.1-10 ng of the ligation product into a chemically competent or electrocompetent cell line that is compatible with the vector.

References

1. Engler, M. J., and Richardson, C. C. (1982) DNA ligases. In *The Enzymes*, Vol. XV (Ed. P. D. Boyer) Academic Press, New York, 3-29.
2. Weiss, B., Thompson, A., and Richardson, C. C. (1968) Enzymatic breakage and joining of deoxyribonucleic acid, VII. Properties of the enzyme-adenylate intermediate in the polynucleotide ligase reaction. *J. Biol. Chem.* 243, 4556-4563.

NxGen® T4 DNA Ligase High Concentration Rapid Kit

Notice of Limited Label License, Copyright, Patents, Warranties, Disclaimers and Trademarks

Lucigen's products are sold for research use only and are not to be used in humans or for medical diagnostics. Lucigen's liability with respect to any T4 DNA Ligase product is limited to the replacement of the product. No other warranties of any kind, expressed or implied, including, without limitation, any implied fitness for any particular use, are provided by Lucigen. Lucigen is not liable for any direct, indirect, incidental or consequential damages arising out of or in connection with the use or inability to use any of its T4 DNA Ligase products.

Some applications in which Lucigen enzymes can be used may be covered by other patents issued and applicable in the United States and certain other countries. Because purchase of this product does not include a license to perform any patented application, users of this product may be required to obtain a patent license depending upon the particular application in which the product is used. It is the sole responsibility of the buyer to ensure that use of the product does not infringe the patent rights of third parties.

Limited Label License

The purchase price of this product includes limited, nontransferable rights to use only the purchased amount of the product. Lucigen Corporation reserves all other rights; in particular, the purchaser of this product may not transfer or otherwise sell this product or its components or derivatives to a third party, and no rights are conveyed to the purchaser to use the product or its components or derivatives for commercial purposes. "Commercial purposes" includes any activity for which a party receives consideration and may include, but is not limited to, (1) use of the product or its components or derivatives in manufacturing, (2) use of the product or its components or derivatives for diagnostic purposes, (3) use of this product or materials made therefrom to provide a service, information, or data (e.g., DNA sequence) to a third party in return for a fee or other consideration, or (4) resale of the product or its components or derivatives, whether or not such product or its components or derivatives are resold for use in research. Academic, Not-for-Profit, and For-Profit institutions must obtain a separate license from Lucigen Corporation to use this product for any purpose other than those permitted above. It is the sole responsibility of the buyer to ensure that use of the product does not infringe that patent rights of third parties.

If the purchaser is not willing to accept these use limitations, Lucigen Corporation is willing to accept return of the product for a full refund. For information on obtaining a license, contact Lucigen Corporation, 2905 Parmenter St., Middleton, WI 53562. Email: Lucigen@lucigen.com. Phone: 608-831-9011. Fax 608-831-9012.