

  
**Lucigen**<sup>®</sup>

Simplifying Genomics


**epicentre**<sup>®</sup>

## RNA 5' Polyphosphatase


Cat. No. RP8092H

1. Introduction

RNA 5' Polyphosphatase\* is a Mg<sup>2+</sup>-independent, 19,000-kDa phosphohydrolase. The enzyme sequentially removes the γ and β phosphates from 5'-triphosphorylated RNA (primary RNA transcripts) and 5'-diphosphorylated RNA.



RNA 5' Polyphosphatase has no activity on RNA with a 5' cap (e.g., 5' m<sup>7</sup>GpppN-----OH 3'), a 5'-monophosphorylated end (5' pN-----OH 3'), or a 5'-hydroxyl end (5' HO-----OH 3'). However, NTPs and dNTPs are used as substrates yielding the corresponding NMPs and dNMPs + inorganic phosphate.



RNA 5' Polyphosphatase is available in a 200 Unit size at a concentration of 20 U/μL. A 10X Reaction buffer is provided with the enzyme.

2. Product Designations and Kit Components

Product	Kit Size	Catalog Number	Reagent Description	Part Numbers	Volume
RNA 5' Polyphosphatase	200 Units	RP8092H	RNA 5' Polyphosphatase (20 Unit/μL)	E0113-20D1	20 μL
			RNA 5' Polyphosphatase 10X Reaction Buffer	SS000807-D1	100 μL

3. Product Specifications

**Storage:** Store only at -20°C in a freezer without a defrost cycle.

**Storage Buffer:** RNA 5' Polyphosphatase is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 1 mM dithiothreitol, 0.1 mM EDTA and 0.1% Triton® X-100.

**Unit Definition:** 1 Unit of RNA 5' Polyphosphatase results in the release of 1 nmol of inorganic phosphate from ATP in 1 hour at 37°C under standard assay conditions.

**RNA 5' Polyphosphatase 10X Reaction Buffer:** 0.5 M HEPES-KOH (pH 7.5), 1 M NaCl, 10 mM EDTA, 1% β-mercaptoethanol, and 0.1% Triton X-100.

**Contaminating Activity Assays:** RNA 5' Polyphosphatase is free of detectable exo- and endonuclease and RNase activities.

**Enzyme Inhibitors:** RNA 5' Polyphosphatase activity is inhibited approximately 50% by 100  $\mu$ M inorganic phosphate, using ATP as a substrate.

#### 4. Applications

- Conversion of 5'-triphosphorylated RNA to 5'-monophosphorylated RNA for use in 5'-T4 RNA Ligase-mediated RNA "tagging" strategies.
- Analysis of 5'-end structure of RNA.
- Preparation of substrate RNA molecules for subsequent degradation using Terminator™ 5'-Phosphate-Dependent Exonuclease (Lucigen).

#### 5. RNA Preparation

1. The RNA should be dissolved in RNase-free Water or RNase-free TE Buffer (10 mM Tris-HCl [pH 7.5], 1 mM EDTA). Total RNA or fractionated RNA preparations can be used.
2. RNA 5' Polyphosphatase activity is inhibited approximately 50% by 100  $\mu$ M inorganic phosphate, using ATP as a substrate. Care should be taken to purify away any residual reaction components so as not to carry over inorganic phosphate into the RNA 5' Polyphosphatase reaction.
3. rNTPs, rNDPs, dNTPs, and dNDPs are all substrates for RNA 5' Polyphosphatase. Therefore, it is strongly recommended that any (deoxy)nucleotide tri- or di- phosphates be removed from the RNA sample before use.

#### RNA 5' Polyphosphatase Standard Protocol

1. Thaw and thoroughly mix the RNA 5' Polyphosphatase 10X Reaction Buffer prior to use.
2. Combine the following components:
 

x	$\mu$ L	RNase-Free water
2	$\mu$ L	RNA 5' Polyphosphatase 10X Reaction Buffer
0.5	$\mu$ L	RiboGuard RNase Inhibitor (optional)
y	$\mu$ L	RNA sample (up to 5 $\mu$ g) (see RNA Preparation)
1	$\mu$ L	RNA 5' Polyphosphatase (20 Units)
<hr/>		
20	$\mu$ L	Total reaction volume
3. Gently but thoroughly mix the reaction.
4. Incubate at 37°C for 30 minutes.
5. Purify the treated RNA by a method appropriate to the downstream application.

*\*RNA 5' Polyphosphatase is covered by patent applications of Epicentre® (an Illumina® company). Purchase of this product is accompanied by a limited, nonexclusive license for the purchaser to use the purchased product solely for lifescience research. Contact Epicentre® (an Illumina® company) for information on licenses for uses in diagnostics or other fields.*

*Epicentre is a trademark of Illumina, Inc. and/or its affiliate(s) in the U.S. and other countries, and is used under license.*

*Terminator is a trademark of Lucigen.*

*Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.*