

sbeadex Lightning buffer recommendations for plant samples

| | Species | Tissue type | Lysis buffer | | | Optimal sbeadex Lightning Core Kit | | | Additional reagents required | |
|---------------|------------|--------------------|---|------------------------------------|--------------------|------------------------------------|--------------------------------|------------------------------|----------------------------------|-----------------------------------|
| | | | Recommended for highest yield | Recommended for highest purity | Non-compatible | Core Kit A (Binding buffer LP) | Core Kit B (Binding buffer LU) | Non-compatible | Protease K Solution ² | Debris capture beads ² |
| Plant samples | Apple | Leaf | LI > H | PN | UR | ✓ | | n/a | ✓ | ✓ |
| | Banana | Leaf | PN > PVP > LI > H | PN | n/a | ✓ | | n/a | ✓ | ✓ |
| | Barley | Seed | PN > BL | PN > UR | n/a | ✓ | ✓ | n/a | ✓ | ✓ |
| | Blackberry | Leaf | UR | UR | PN, PVP, BL, H, LI | | ✓ | n/a | ✓ | ✓ |
| | | Leaf | UR > PN > LI | UR > PN | BL | ✓ | | n/a | ✓ | ✓ |
| | Canola | Seed | LI > UR | LI > UR | n/a | ✓ | | n/a | ✓ | ✓ |
| | Carrot | Leaf | UR > PN > PVP | UR > PN > PVP | BL, H, LI | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ |
| | | Leaf | BL > LBH | PN > UR | n/a | | ✓ | n/a | ✓ | ✓ |
| | Corn | Seed | UR > LI > H > BL | PN | n/a | ✓ | | n/a | ✓ | ✓ |
| | | Leaf | LI > PVP | LI > VP | n/a | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ |
| | Cucumber | Seed | LI | H | PN, PVP, BL, UR | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ |
| | Currant | Leaf | PVP > PN | PVP > PN | BL, H, UR, LI | ✓ | | n/a | ✓ | ✓ |
| | Hop | Leaf | LI > UR | UR > PN > PVP > LI | BL, H | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ |
| | Lentil | Seed | LI > PVP | LI > PVP | BL, H | ✓ | | n/a | ✓ | ✓ |
| | | Leaf | LI > H > UR > BL > PN > PVP | H > UR > PN > PVP | n/a | ✓ | | n/a | ✓ | ✓ |
| | Oat | Seed | LI > H > BL | LI > BL > H | n/a | ✓ | ✓ | n/a | ✓ | ✓ |
| | Pea | Seed | LI > PVP | LI > PVP > PN | BL, H, UR | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ |
| | Potato | Leaf | UR > PN | UR > PN | BL | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ |
| | | Leaf | UR > BL > LI | UR > PN > PVP > LI | n/a | ✓ | | n/a | ✓ | ✓ |
| | Rice | Seed | BL > LI > H | LI > H > UR | n/a | ✓ | | n/a | ✓ | ✓ |
| | Sorghum | Seed | BL > UR > LI > PN | PN > UR > LI | n/a | ✓ | ✓ | n/a | ✓ | ✓ |
| | | Leaf | LI > H > BL > UR | PN > UR | n/a | ✓ | | n/a | ✓ | ✓ |
| | Soy | Seed | LI > H > BL | LI | PN, PVP | | ✓ | n/a | ✓ | ✓ |
| | Sugar beet | Leaf | UR > LI > BL | UR > LI > H > PN | n/a | ✓ | | n/a | ✓ | ✓ |
| | Sunflower | Leaf | UR > PN | UR > PN | PVP, BL, H, LI | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ |
| | | Leaf | UR > LI | UR > LI > H > PN | n/a | ✓ | | n/a | ✓ | ✓ |
| Tomato | Seed | UR > PVP | UR > PVP | BL, H, LI | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ | |
| | Leaf | LI > UR | LI > UR | n/a | ✓ | | LU (Core Kit B) ¹ | ✓ | ✓ | |
| Wheat | Seed | UR > LI | PN | n/a | ✓ | | n/a | ✓ | ✓ | |

Table 1. Lysis buffer performance for plant tissue samples.

LGC Biosearch Technologies have tested a range of plant species with [sbeadex Lightning chemistry](#) to determine the optimal lysis buffer for DNA yield and for DNA purity. Lysis buffers tested with each plant sample type were PN, PVP, BL, H, UR and LI. The table details the optimal lysis buffer for each sample type (in bold). All alternative lysis buffers are then listed for either yield or purity in descending order of suitability, and non-compatible lysis buffers are listed separately. The most appropriate sbeadex Lightning core kit for each sample type is also detailed. Note that buffers and core kits listed as incompatible for a particular sample type are based on our in-house testing; it is possible that with sample variation (e.g. younger leaves, alternative input volume), different results may be obtained.

Biosearch Technologies have provided this information as guidance but strongly recommend that you use the [sbeadex Lightning Nucleic Acid Purification Starter Kit](#) to ensure the optimal nucleic acid purification conditions for your specific sample type and downstream application.

For more detailed information regarding this data, or for any other technical queries regarding sbeadex Lightning, please contact our technical support team: techsupport@lqcgroupp.com.

¹for some lysis buffers

²2 µL Protease K solution (20 mg/mL) per 1 mL of lysis buffer. 40 µL Debris capture beads per 1 mL of lysis buffer.

sbeadex Lightning buffer recommendations for livestock samples

| | Species | Tissue type | Lysis buffer | | | Optimal sbeadex Lightning Core Kit | | | Additional reagents required | |
|-------------------|---------|--------------------------------------|-----------------------|--------------------|------------------|------------------------------------|--------------------------------|----------------|----------------------------------|-----------------------------------|
| | | | Recommended for yield | Optimal for purity | Non-compatible | Core Kit A (Binding buffer LP) | Core Kit B (Binding buffer LU) | Non-compatible | Protease K Solution ⁵ | Debris capture beads ⁵ |
| Livestock samples | Bovine | Muscle tissue | LI | LI | n/a | ✓ | | n/a | ✓ | |
| | | Hair | H | LI | n/a | ✓ | | n/a | ✓ | |
| | | Ear punch (in preservative solution) | PVP ¹ | PVP ³ | n/a | ✓ | | n/a | ✓ | |
| | Chicken | Muscle tissue | LI | LI | n/a | ✓ | | n/a | ✓ | |
| | Fish | Tissue | LI | LI | n/a ⁴ | | ✓ | n/a | ✓ | |

Table 2. Lysis buffer performance for livestock tissue samples.

Biosearch Technologies have tested a range of livestock species with [sbeadex Lightning chemistry](#) to determine the optimal lysis buffer for DNA yield and for DNA purity. Not all sbeadex Lightning lysis buffers were tested with all sample types due to known incompatibilities with sample preservation solutions. The table details the optimal lysis buffer for each sample type. Incompatibility of lysis buffers with livestock sample types is highly dependent on sample type; the buffers listed as optimal for yield and for purity are recommended for these sample types. The most appropriate sbeadex Lightning core kit for each sample type is also detailed in the table.

Biosearch Technologies have provided this information as guidance but strongly recommend that you use the [sbeadex Lightning Nucleic Acid Purification Starter Kit](#) to ensure the optimal nucleic acid purification conditions for your specific sample type and downstream application. If working with samples in Allflex buffer, or using Lysis buffer PVP, an additional wash step with Wash Buffer TN2 or PN2 is strongly recommended.

For more detailed information regarding this data, or for any other technical queries regarding sbeadex Lightning, please contact our technical support team: techsupport@lgcgroup.com.

³Wash with Wash buffer TN2 instead of water

⁴Lysis buffer BL not tested for this sample type

⁵10 µL Protease K solution per sample. 40 µL Debris capture beads per 1 mL of lysis buffer.