

Manual

Exonuclease I, *E. coli*

For Research Use Only. Not for use in diagnostic procedures.

Exonuclease I, *E. coli* is part of the Epicentre™ product line, known for its unique genomics kits, enzymes, and reagents which offer high quality and reliable performance.

Manual

Exonuclease I, *E. coli*

Contents

1. Introduction	3
2. Product designations and kit components	3
3. Product specifications	3
4. Applications	3
5. References	4
6. Further support	4

Manual

Exonuclease I, *E. coli*

1. Introduction

Exonuclease I (Exo I), the product of the *sbcB* gene of *E. coli*, is an exodeoxyribonuclease that hydrolyses single-stranded DNA (ssDNA) stepwise in a 3'→5' direction.¹⁻³ Hydrolysis generates deoxyribonucleoside 5' monophosphates and a terminal dinucleotide diphosphate.¹ The enzyme requires magnesium (optimal Mg²⁺ concentration is 10 mM) and the presence of a free 3'-hydroxyl terminus.¹ Exo I is active under a wide variety of buffer conditions, allowing addition of the enzyme directly into most reaction mixes. Heat inactivation results from incubation at 80 °C for 15 minutes.

2. Product designations and kit components

Product	Kit size	Catalog number	Reagent description	Part number	Volume
Exonuclease I, <i>E. coli</i>	20,000 Units	X40520K	Exonuclease I (20 Units/μL)	E0027-20D1	1 mL

3. Product specifications

Storage: Store only at -20 °C in a freezer without a defrost cycle.

Storage buffer: Exo I is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 1.0 mM dithiothreitol, 0.1 mM EDTA and 0.1% Triton® X-100 (Rohm & Haas).

Unit definition: One unit of Exo I results in the acid-solubilisation of 10 nmol of nucleotides from calf thymus DNA in 30 minutes at 37 °C.

Quality control: Exo I is function-tested in a reaction containing 33 mM Tris-acetate (pH 7.5), 66 mM potassium acetate, 10 mM magnesium acetate, 0.5 mM dithiothreitol, 10 μg of denatured calf thymus DNA and varying amounts of Exo I.

Contaminating activity assays: Exo I is free of detectable RNase, endonuclease and double-stranded exonuclease activities.

4. Applications

Removal of residual ssDNA and oligonucleotides from reaction mixes. Linear ssDNA and oligonucleotides can be selectively degraded from heterogeneous mixtures of nucleic acids in reaction mixes.

Clean-up of PCR. Unused amplification primers can be removed after PCR using Exo I.

Removal of ssDNA from nucleic acid mixtures. Linear ssDNA can be selectively degraded from heterogeneous mixtures of nucleic acids with Exo I.

Assay for regions of ssDNA.^{4,5} Use Exo I to assay for the presence of ssDNA containing a free 3'-hydroxyl end. This technique was used to characterise the endonuclease and helicase activities of purified *recBC* protein on circular fd phage DNA and duplex phage T7 DNA respectively.

Manual

Exonuclease I, *E. coli*

5. References

1. Lehman, I.R. and Nussbaum, A.L. (1964) *J. Biol. Chem.* **239**, 2628.
2. Kusher, S.R. et al., (1971) *Proc. Natl. Acad. Sci. USA* **68**, 824.
3. Kusher, S.R. et al., (1972) *Proc. Natl. Acad. Sci. USA* **69**, 1366.
4. Goldmark, P.J. and Linn, S. (1972) *J. Biol. Chem.* **247**, 1849.
5. Rosamond, J. et al., (1979) *J. Biol. Chem.* **254**, 8646.

6. Further support

If you require any further support, please do not hesitate to contact our Technical Support Team:

techsupport@lgcgroup.com



   @LGCBiosearch | biosearchtech.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2021. All rights reserved. GEN/888/EK/0221

For Research Use Only. Not for use in diagnostic procedures.



BIOSEARCH™
TECHNOLOGIES

GENOMIC ANALYSIS BY LGC