

# Manual

## Endonuclease VIII Kit

---

For Research Use Only. Not for use in diagnostic procedures.

# Manual

Endonuclease VIII Kit

---

## Contents

<a href="#">1. Introduction .....</a>	<a href="#">3</a>
<a href="#">2. Product designations and kit components .....</a>	<a href="#">3</a>
<a href="#">3. Product specifications.....</a>	<a href="#">3</a>
<a href="#">4. References.....</a>	<a href="#">4</a>

# Manual

## Endonuclease VIII Kit

### 1. Introduction

*E. coli* Endonuclease VIII is a bifunctional DNA glycosylase with DNA N-glycosylase and AP lyase activities. The N-glycosylase activity releases damaged pyrimidines from double-stranded DNA generating an apurinic (AP) site. The AP lyase activity cleaves the DNA phosphodiester backbone at AP sites via  $\beta$ - and  $\delta$ -elimination, creating a one nucleotide gap with 5' and 3' phosphate termini. Damaged bases recognised and removed by Endonuclease VIII include thymine glycol, uracil glycol, urea, 5,6-dihydroxythymine, 5-hydroxy-5-methylhydantoin, 5-hydroxy-6-hydrothymine, 5,6-dihydrouracil and alloxan. Although Endonuclease VIII is like Endonuclease III, Endonuclease VIII has both  $\beta$  and  $\delta$  lyase activity while Endonuclease III has only  $\beta$  lyase activity..

### 2. Product designations and kit components

Product	Kit size	Catalogue number	Reagent description	Part number	Volume
Endonuclease VIII Kit	1,000 units	30270	Endonuclease VIII (10 U/ $\mu$ L)	F836129-1	100 $\mu$ L
			Endonuclease VIII Reaction Buffer, 10X	F836075-1	1 mL
	10,000 units	30271	Endonuclease VIII (10 U/ $\mu$ L)	F836129-2	1 mL
			Endonuclease VIII Reaction Buffer, 10X	F836075-2	5 mL

### 3. Product specifications

**Storage:** Store at -20 °C in a freezer without a defrost cycle.

**Storage buffer:** Endonuclease VIII is supplied in a 50% glycerol solution containing 10 mM Tris-HCl (pH 8.0), 250 mM NaCl, 0.1 mM EDTA.

**Endonuclease VIII Reaction Buffer, 10X:** 100 mM Tris-HCl (pH 8.0), 750 mM NaCl, 10 mM EDTA.

**Unit definition:** One unit is the amount of enzyme required to cleave 1 pmol of 72 base pair oligonucleotide duplex containing a single abasic\* site on each strand in 1 hour at 37 °C in 25  $\mu$ L in 1x Endonuclease VIII reaction buffer containing 10 pmol of oligonucleotide duplex.

\*Abasic site is created by including Uracil N-Glycosylase (UNG) concurrently with the Endonuclease VIII.

**Contaminating activity assays:** Endonuclease VIII is free of detectable Endonuclease, Exonuclease and RNase activities.

**Purity:** Endonuclease VIII is assessed as  $\geq$ 95% pure by SDS-Polyacrylamide gel electrophoresis.

**Heat inactivation:** Endonuclease VIII can be inactivated by heating at 75 °C for 10 minutes.

### 4. References

# Manual

## Endonuclease VIII Kit

---

Kropachev, K. Y. *et al.* (2006) *Biochemistry*, 45, 12039-12049. The catalytic mechanism of *E. coli* endonuclease VIII: Roles of the intercalation loop and the zinc finger.

Pflaum, M. *et al.* (1998) *Free Rad. Res*, 29, 585-594. Alkaline elution. DNA oxidation products determined with repair endonucleases in mammalian cells: Types, basal levels and influence of cell proliferation.

Jiang, D. *et al.* (1997) *J. Biol. Chem.*, 272, 32230-32239. Characterization of *Escherichia coli* Endonuclease VIII.

Hartwig, *et al.* (1996). *Toxicol. Lett*, 88, 85-90. Alkaline unwinding. Sensitive analysis of oxidative DNA damage in mammalian cells.

Collins, A. *et al.* (1996). *Environ. Health Perspect*, 104, 465-469. Comet assay. Direct enzymatic detection of endogenous oxidative base damage in human lymphocyte DNA.

Demple, B, and Harrison, L. (1994) *Annu. Rev. Biochem*, 63, 915-48 (see page 923). Repair of Oxidative Damage to DNA: Enzymology and Biology.

Hatahet, Z. *et al.* (1994) *J. Biol. Chem.*, 269, 18814-18820. New substrates for old enzymes.

Dizdaroglu, M. *et al.* (1993) *Biochemistry* 32, 12105-12111. Substrate specificity of the *Escherichia coli* endonuclease III: Excision of thymine- and cytosine-derived lesions in DNA produced by radiation-generated free radicals.

Collins, A. *et al.* (1993). *Carcinogenesis*, 14, 1733-1735. Comet assay. A simple technique for quantitation of low levels of DNA damage in individual cells.

Singh, N. *et al.* (1988) *Exp. Cell. Res*, 175, 184-191. Comet assay. Oxidative damage to DNA: Do we have a reliable biomarker?



   @LGCBiosearch | biosearchtech.com

*All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2023. All rights reserved. GEN/1043/EK/0123*

For Research Use Only. Not for use in diagnostic procedures.



**BIOSEARCH™**  
**TECHNOLOGIES**

GENOMIC ANALYSIS BY LGC